



abaxis healthcare



abachem 300

MORE VALUABLE THAN USUAL

Premium Auto Biochemistry Analyzer

Swipe card reader for florentine ASO,RF,CRP & HbA1C Tests


5" Graphic LCD with colour touch panel

Real time online kinetic graph

Isolated circuit for Lamp protection

Monochromatic and bichromatic measurement

Multi Standard Calibration



abachem 300

Clinical Chemistry Analyzer

TECHNICAL SPECIFICATIONS

Human Machine Interface	: Touch Panel/ Keypad /EXTERNAL KEY BOARD
Linear Measurement Range	: 0.000 to 3.000 Absorbance unit (A)
Photometric Accuracy	: $\pm 2\%$ or 0.07 whichever is higher, from 0 to 1.5 A $\pm 3\%$ from 1.5 A to 3.0 A
Drift	: <0.007 A/hr.
Photometric range	: 2.2A
Optical Measurement	: Photodiode

FILTERS

Type of Filter	: High quality narrow band interference
Wave Length	: 340, 405, 510, 546, 578 and 630 nm & Two Optional
Half bandwidth	: 10nm \pm 2nm
Filter Selection	: Automatic by Stepper Motor

FLOW CELL / CUVETTE

Sipping Volume	: 300 to 1000 μ l
Flow Cell Volume	: 18 μ l
Sipping mode	: Automatic by specially designed Peristaltic Pump
Cuvette Volume	: 500 μ l
Temperature of Cuvette / flow cell Block	
Method	: By Peltier effect
Temperature	: 25° C, 30° C and 37° C

Light Source	: Tungsten Halogen
Warm up Time	: 90 Sec.
Display	: 5" Graphics LCD, Negative Blue, STN
Printer	: Built in thermal printer, 20 Column

Memory	: 64 KB - Non volatile RAM with Battery back up
Storage Capacity	: 250 Open tests with 30 QC results/test (Normal and abnormal controls) and 2500 Patient results with patient ID

ANALYSIS MODE	: Absorbance	End point
	: Fixed time	Kinetic
	: Differential	Ratio,
	: Coagulation (PT, INR)	e-c mode

Concentration calculation	: By factor or by standard or by electronic calibration
RS232 serial port	: 9600 baud, 1 start, 8 data, 1 stop, no parity bits

POWER

Wattage	: 50 Watts
Voltage	: 115 - 230 Volts \pm 10 %, 50 - 60 Hz\

Operating Position	: On horizontal, flat, rigid and vibration free surface Operating Conditions
--------------------	------------------------------------------------------------------------------

Temperature	: from + 18° C to 35° C
Relative Humidity	: up to 85 %

STORAGE CONDITION

Temperature	: from - 10° C to + 60° C
Relative Humidity	: up to 85 %
Enclosure	: ABS Fire retardant
Size (cm)	: 30 X 38 X 13.5 (l x b x h)
Weight	: 5.5 Kgs. (Approx)

Ordering Information :

Product Code : F1010-0162

Abaxis Healthcare India Pvt Ltd
#134, Gandhi Nagar main, Cheror Post Thrissur
Kerala - 680 008, Phone: 0487 2323 537, +91 9961 855 000
e-mail: sales@abaxishealthcare.com, www.abaxishealthcare.com



abaxis healthcare